



Barostim
Outsmart the heart

Surgical Procedures

Caution: Federal law restricts this device to sale by or on the order of a physician

ABOUT THIS DOCUMENT

This document is a portion of the Instructions for Use (IFU) for the Barostim™ System. The full IFU consists of:

System Overview	900133-001
Surgical Procedures	900133-002
Programming	900133-003
Magnetic Resonance Imaging (MRI)	900133-004
Patient Instructions	900133-005

IFU documents are available at www.cvr.com/ifu.

Table of Contents

About this document.....	2
1 Implant Preparation	4
WARNINGS and PRECAUTIONS	5
Before Implantation	5
Materials Recommended for Implantation	5
Materials Needed	6
Surgical Preparation.....	6
2 Implantation Procedure	9
Clinical Update.....	10
Anesthesia	11
Lead Placement.....	13
IPG Placement.....	13
Connect the Lead	14
Mapping & Fixating the Lead.....	16
Final positioning of the lead and Fixation	18
3 IPG Replacement.....	21
Recommendations	22
Materials Needed	22
Antibiotic Coverage.....	22
Explantation of Depleted IPG	22
IPG Placement.....	23
4 Explantation Procedures.....	24
Recommendations	25
Materials Recommended	25
Antibiotic Coverage – Not Infected	25
Antibiotic Coverage – Infected	25
IPG Explantation	25
CSL Explantation	26

1

Implant Preparation

WARNINGS AND PRECAUTIONS

For a complete listing of Warnings and Precautions for the system including implant procedures please reference System Overview 900133-001 (available at www.cvr.com/ifu).

BEFORE IMPLANTATION

Formal preoperative duplex ultrasonography should.

- Confirm the absence of complex arterial anatomy, such as carotid kinks, loops and coils, which would compromise the implant procedure.
- Verify absence of any stenosis producing a greater than 50% reduction in diameter of the carotid arteries.
- Verify absence of any ulcerative plaques.
- Verify the level of the carotid bifurcation is easily accessible from standard cervical incisions.
- Determine if there are any anatomic variants present which might suggest that additional imaging would be helpful for treatment planning.
- Other variants or conditions that would eliminate the patient as a surgical candidate.

Ensure that a backup Implantable Pulse Generator (IPG) and Carotid Sinus Lead (CSL or Lead) are available in case the sterility of the first IPG is compromised, or the device is damaged during surgery.

Ensure that a backup Programmer System is available in case the primary system is damaged or becomes non-operational.

It is recommended that physicians verify compatibility with the implanted device(s) at the time of Barostim System implant. Ensure that proper instrumentation for monitoring the behavior of the currently implanted device is present.

MATERIALS RECOMMENDED FOR IMPLANTATION

- A table or stand outside the sterile operating field to hold the Programmer System.
- Blood pressure monitoring equipment (such as an arterial line) for assessment of blood pressure changes during therapy testing.